



by MAVIG

MCatalogue & MPortfolio 2025/1

Index

M

General Information



RA431

All-round Protection – Costume (Two-piece set: vest & skirt)

Article numbers

2



RA432

All-round Protection – Coat

Article numbers

3



RA460

Front Protection – Apron

Article numbers

4

10+11



RA410

Head Protection – Cap

Article numbers

5

12



RA414

Thyroid Protection

Article numbers

6

13



RA437

Shoulder & Chest Protection – Bolero

Article numbers

7

14



AW

Storage Systems

Article numbers

7

14

7

14

8

15

Minimalistic. MAVIG. Munich.

Our budget solution.

„M by MAVIG“ is the new, reduced product line of MAVIG X-ray protective clothing for users – Made in Germany.

With our „M“-Selection (3 models, 4 lengths, 4 colours, 5 sizes and accessories) you can quickly find effective radiation protection for every need.

„M by MAVIG“ complies with the standards and certifications.

M by MAVIG: Where safety meets cost efficiency.

MCert

Standard-compliant

that means:

- IEC 61331-3:2014 / DIN EN 61331-3:2016
- CE marking according to PPE regulation
- MDR (EU 2017/745)
- Protection over 60% of the largest circumference of the user's body
- and much more

MPure

X-ray protective material

- Light lead-containing material
- Best possible combination of weight and radiation protection
- Radiation protection in the entire X-ray tube voltage range from **50 - 150 kV**
- Quality according to standards (IEC 61331-1:2014 / DIN EN 61331-1:2016)

MFree

X-ray protective material

- Extremely light lead-free material
- Optimised for long wearing times
- Radiation protection in the entire X-ray tube voltage range from **50 - 110 kV**
- Quality according to standards (IEC 61331-1:2014 / DIN EN 61331-1:2016)

MTex

Outer material

- Durable, easy to clean nylon fabric
- Special PU coating
- OEKOTEX 100 Standard

Colours

- 51 Inky Blue
- 52 New Red
- 53 Green 17
- 54 Modern Grey

RA431 All-round protection Costume (vest & skirt)



FRONT



BACK



SIDE

Two-piece set – Vest & Skirt

All-round protection with ergonomic comfort thanks to the division into vest and skirt

- Full lead equivalent in both front wings, no full-surface overlapping necessary
- Cut to fit the body
- Especially suitable for long wearing times
- Weight relief for the spine
- Shoulder pads as standard
- Type examination and certification according to (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

Stretch belt for skirt, optional



Protective material/Lead equivalent

MPure – Light lead (50 - 150kV)
or MFree – Lead free (50 - 110kV)

Front panel: 0.25 mm Pb,
0.35 mm Pb
or 0.50 mm Pb

Generally in the back: 0.25 mm Pb

Outer material/Colours

MTE – Nylon fabric

Our colours:

	Inky Blue
	New Red
	Green 17
	Modern Grey

Options

- Additional stretch belt (skirt)
- Sewn-on shoulder border with embroidery, one line of text only
- Breast pocket, only with logo embroidery (outer regular insertion pocket)

Sizing according to standard

The person's height, as well as either circumference of the chest, waist or hip, depending on which of the three dimensions is the greatest, is used for this purpose.

RA431 Sizes

	XS	S	M	L	XL
Breast (for vest)	≤ 84 cm	85 - 95 cm	96 - 106 cm	107 - 117 cm	118 - 128 cm
Waist/hip (for skirt)	≤ 86 cm	87 - 97 cm	98 - 107 cm	108 - 118 cm	119 - 129 cm

Lengths (skirt)

for a body height of:	65 cm	70 cm	75 cm
	150 - 170 cm	170 - 180 cm	180 - 190 cm

RA432 All-round protection Coat



FRONT



BACK



SIDE

X-ray protective Coat

All-round protection for quick put-on thanks to the one-piece design

- Full lead equivalent in both front wings, no full-surface overlapping necessary
- Cut to fit the body
- Effective shoulder relief thanks to the elastic, exchangeable stretch belt
- Well suited for long wearing times
- Shoulder pads as standard
- Type examination and certification according to (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



Protective material/Lead equivalent

MPure – Light lead (50 - 150kV)
or MFree – Lead free (50 - 110kV)

Front panel: 0.25 mm Pb,
0.35 mm Pb
or 0.50 mm Pb

Generally in the back: 0.25 mm Pb

Outer material/Colours

MText – Nylon fabric

Our colours:

	Inky Blue
	New Red
	Green 17
	Modern Grey

Options

- Sewn-on shoulder border with embroidery, one line of text only
- Breast pocket, only with logo embroidery
- (outer regular insertion pocket)

Sizing according to standard

The person's height, as well as either circumference of the chest, waist or hip, depending on which of the three dimensions is the greatest, is used for this purpose.

RA432 Sizes

	XS	S	M	L	XL
Breast	≤ 84 cm	85 - 95 cm	96 - 106 cm	107 - 117 cm	118 - 128 cm
Waist/hip	≤ 98 cm	99 - 108 cm	109 - 119 cm	120 - 130 cm	131 - 141 cm

Lengths (coat)

for a body height of:	90 cm	100 cm	110 cm	120 cm
	150 - 160 cm	160 - 170 cm	170 - 180 cm	175 - 190 cm

RA460 Front protection Apron



FRONT



BACK



SIDE

X-ray protective Apron

Lightweight front protection for applications with radiation exposure from the front

- Snap buckle fastener at the front for easy fastening
- Breathable mesh in the back wings
- Shoulder pads as standard
- Type examination and certification according to (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016



Protective material/Lead equivalent

MPure – Light lead (50 - 150kV)
or MFree – Lead free (50 - 110kV)

Front panel: 0.25 mm Pb,
0.35 mm Pb
or 0.50 mm Pb
Backside: Mesh fabric

Outer material/Colours

MTE – Nylon fabric

Our colours:

	Inky Blue
	New Red
	Green 17
	Modern Grey

Options

- Sewn-on shoulder border with embroidery, one line of text only
- Breast pocket, only with logo embroidery (outer regular insertion pocket)

Sizing according to standard

The person's height, as well as either circumference of the chest, waist or hip, depending on which of the three dimensions is the greatest, is used for this purpose.


RA460 Sizes

	XS	S	M	L	XL
Breast	≤ 84 cm	85 - 95 cm	96 - 106 cm	107 - 117 cm	118 - 128 cm
Waist/hip	≤ 91 cm	92 - 102 cm	103 - 112 cm	113 - 123 cm	124 - 134 cm

Lengths (apron)

	90 cm	100 cm	110 cm	120 cm
for a body height of:	150 - 160 cm	160 - 170 cm	170 - 180 cm	175 - 190 cm


RA410 Head protection



Radiation protection cap to protect the skull

- Breathable mesh top
- Universal size: suitable for all head circumferences from 54 - 63 cm
- Exclusively available in Modern Grey
- Type examination and certification according to (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

RA414 Thyroid protection



Thyroid sternum protection to protect the particularly radiation-sensitive thyroid glands

- Easy to put on thanks to hook&burr fastener
- Universal size: suitable for all neck circumferences from 32 - 48 cm
- Fits the whole „M“ products
- Exclusively available in Modern Grey
- Type examination and certification according to (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

RA437 Shoulder & chest protection – Bolero



Radiation protection for the shoulder joint and breast tissue without restricting freedom of movement

- Separate radiation protection compatible with all models
- Closure with hook & burr fastener in the front area
- Sizes matching the X-ray protective apron
- Exclusively available in Modern Grey
- Type examination and certification according to (EU) 2016/425, (EU) 2017/745, IEC 61331-1:2014 / DIN EN 61331-1:2016 and IEC 61331-3:2014 / DIN EN 61331-3:2016

AW106 Special hanger



AW106

Clothes hanger made from heavy-duty material (100% recycled) for the storage of individual aprons, coats or costumes

- 100% recycling material
- for X-ray protective clothing (1 pc.)
- Dimensions: 90 x 520 mm (H x W)

AW415 Wall system



AW41502



AW41503




AW41506

Proper storage of two, three or six items of clothing for wall mounting

- Storage on special hangers (AW106)
- Convenient access
- Proper storage
- Low space requirement
- Easy assembly and install

AW425 Mobile system



AW42506



AW42512

Mobile system for the proper storage of six or twelve items of clothing

- Storage on special hangers (AW106)
- High mobility and flexibility
- Proper storage

Further storage systems can be found on our website:

www.mavig.com

Article numbers.

*Article numbers for our „M by MAVIG“ product line
can be found on the following pages*



M by MAVIG : Where safety meets cost efficiency.



Costume

	Front (Pb)	Lenght (Skirt)	Article number* (Start)	(Colour code)	MPure (Inner material)	MFree (Inner material)
X-Small (XS)	0.25 mm Pb	65 cm	RA431K22Z65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431K22Z70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431K22Z75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431K32Z65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431K32Z70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431K32Z75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431K52Z65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431K52Z70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431K52Z75	+CC	+55	+56
Small (S)	0.25 mm Pb	65 cm	RA431K22Z65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431K22Z70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431K22Z75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431K32S65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431K32S70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431K32S75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431K52S65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431K52S70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431K52S75	+CC	+55	+56
Medium (M)	0.25 mm Pb	65 cm	RA431K22M65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431K22M70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431K22M75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431K32M65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431K32M70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431K32M75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431K52M65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431K52M70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431K52M75	+CC	+55	+56
Large (L)	0.25 mm Pb	65 cm	RA431K22L65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431K22L70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431K22L75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431K32L65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431K32L70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431K32L75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431K52L65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431K52L70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431K52L75	+CC	+55	+56
X-Large (XL)	0.25 mm Pb	65 cm	RA431K22X65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431K22X70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431K22X75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431K32X65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431K32X70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431K32X75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431K52X65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431K52X70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431K52X75	+CC	+55	+56

Colour codes



Inky Blue	51
New Red	52
Green 17	53
Modern Grey	54

Optional

Additional stretch belt (skirt) - Size XS - S
 Additional stretch belt (skirt) - Size M - L
 Additional stretch belt (skirt) - Size XL
 Sewn-on shoulder border with embroidery (one line of text)
 Breast pocket (outer regular insertion pocket) without embroidery
 Breast pocket (outer regular insertion pocket) with embroidery
 Breast pocket (outer regular insertion pocket) with logo embroidery

* Example article number: RA431K52M705255 = Costume with 0.50 mm Pb (Front), size Medium, skirt length 70 cm, in New Red and with MPure, MText

Vest

	Front (Pb)	Article number (Start)	(Colour code)	MPure (Inner material)	MFree (Inner material)
XS	0.25 mm Pb	RA431W22Z00	+CC	+55	+56
	0.35 mm Pb	RA431W32Z00	+CC	+55	+56
	0.50 mm Pb	RA431W52Z00	+CC	+55	+56
S	0.25 mm Pb	RA431W22S00	+CC	+55	+56
	0.35 mm Pb	RA431W32S00	+CC	+55	+56
	0.50 mm Pb	RA431W52S00	+CC	+55	+56
M	0.25 mm Pb	RA431W22M00	+CC	+55	+56
	0.35 mm Pb	RA431W32M00	+CC	+55	+56
	0.50 mm Pb	RA431W52M00	+CC	+55	+56
L	0.25 mm Pb	RA431W22L00	+CC	+55	+56
	0.35 mm Pb	RA431W32L00	+CC	+55	+56
	0.50 mm Pb	RA431W52L00	+CC	+55	+56
XL	0.25 mm Pb	RA431W22X00	+CC	+55	+56
	0.35 mm Pb	RA431W32X00	+CC	+55	+56
	0.50 mm Pb	RA431W52X00	+CC	+55	+56

Skirt

	Front (Pb)	Lenght (Skirt)	Article number (Start)	(Colour code)	MPure (Inner material)	MFree (Inner material)
X-Small (XS)	0.25 mm Pb	65 cm	RA431R22Z65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431R22Z70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431R22Z75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431R32Z65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431R32Z70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431R32Z75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431R52Z65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431R52Z70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431R52Z75	+CC	+55	+56
Small (S)	0.25 mm Pb	65 cm	RA431R22S65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431R22S70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431R22S75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431R32S65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431R32S70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431R32S75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431R52S65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431R52S70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431R52S75	+CC	+55	+56
Medium (M)	0.25 mm Pb	65 cm	RA431R22M65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431R22M70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431R22M75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431R32M65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431R32M70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431R32M75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431R52M65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431R52M70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431R52M75	+CC	+55	+56
Large (L)	0.25 mm Pb	65 cm	RA431R22L65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431R22L70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431R22L75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431R32L65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431R32L70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431R32L75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431R52L65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431R52L70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431R52L75	+CC	+55	+56
X-Large (XL)	0.25 mm Pb	65 cm	RA431R22X65	+CC	+55	+56
	0.25 mm Pb	70 cm	RA431R22X70	+CC	+55	+56
	0.25 mm Pb	75 cm	RA431R22X75	+CC	+55	+56
	0.35 mm Pb	65 cm	RA431R32X65	+CC	+55	+56
	0.35 mm Pb	70 cm	RA431R32X70	+CC	+55	+56
	0.35 mm Pb	75 cm	RA431R32X75	+CC	+55	+56
	0.50 mm Pb	65 cm	RA431R52X65	+CC	+55	+56
	0.50 mm Pb	70 cm	RA431R52X70	+CC	+55	+56
	0.50 mm Pb	75 cm	RA431R52X75	+CC	+55	+56

Colour codes

	Inky Blue	51
	New Red	52
	Green 17	53
	Modern Grey	54

Coat

	Front (Pb)	Length (Coat)	Article number* (Start)	(Colour code)	MPure (Inner material)	MFree (Inner material)
X-Small (XS)	0.25 mm Pb	90 cm	RA43222Z90	+CC	+55	+56
	0.25 mm Pb	100 cm	RA43222Z10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA43222Z11	+CC	+55	+56
	0.35 mm Pb	90 cm	RA43232Z90	+CC	+55	+56
	0.35 mm Pb	100 cm	RA43232Z10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA43232Z11	+CC	+55	+56
	0.50 mm Pb	90 cm	RA43252Z90	+CC	+55	+56
	0.50 mm Pb	100 cm	RA43252Z10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA43252Z11	+CC	+55	+56
Small (S)	0.25 mm Pb	90 cm	RA43222S90	+CC	+55	+56
	0.25 mm Pb	100 cm	RA43222S10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA43222S11	+CC	+55	+56
	0.35 mm Pb	90 cm	RA43232S90	+CC	+55	+56
	0.35 mm Pb	100 cm	RA43232S10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA43232S11	+CC	+55	+56
	0.50 mm Pb	90 cm	RA43252S90	+CC	+55	+56
	0.50 mm Pb	100 cm	RA43252S10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA43252S11	+CC	+55	+56
Medium (M)	0.25 mm Pb	100 cm	RA43222M10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA43222M11	+CC	+55	+56
	0.25 mm Pb	120 cm	RA43222M12	+CC	+55	+56
	0.35 mm Pb	100 cm	RA43232M10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA43232M11	+CC	+55	+56
	0.35 mm Pb	120 cm	RA43232M12	+CC	+55	+56
	0.50 mm Pb	100 cm	RA43252M10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA43252M11	+CC	+55	+56
	0.50 mm Pb	120 cm	RA43252M12	+CC	+55	+56
Large (L)	0.25 mm Pb	100 cm	RA43222L10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA43222L11	+CC	+55	+56
	0.25 mm Pb	120 cm	RA43222L12	+CC	+55	+56
	0.35 mm Pb	100 cm	RA43232L10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA43232L11	+CC	+55	+56
	0.35 mm Pb	120 cm	RA43232L12	+CC	+55	+56
	0.50 mm Pb	100 cm	RA43252L10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA43252L11	+CC	+55	+56
	0.50 mm Pb	120 cm	RA43252L12	+CC	+55	+56
X-Large (XL)	0.25 mm Pb	100 cm	RA43222X10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA43222X11	+CC	+55	+56
	0.25 mm Pb	120 cm	RA43222X12	+CC	+55	+56
	0.35 mm Pb	100 cm	RA43232X10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA43232X11	+CC	+55	+56
	0.35 mm Pb	120 cm	RA43232X12	+CC	+55	+56
	0.50 mm Pb	100 cm	RA43252X10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA43252X11	+CC	+55	+56
	0.50 mm Pb	120 cm	RA43252X12	+CC	+55	+56



Colour codes



Inky Blue	51
New Red	52
Green 17	53
Modern Grey	54

Optional

Sewn-on shoulder border with embroidery (one line of text)
 Breast pocket (outer regular insertion pocket) without embroidery
 Breast pocket (outer regular insertion pocket) with embroidery
 Breast pocket (outer regular insertion pocket) with logo embroidery

* Example article number: RA43252M125355 = Coat with 0.50 mm Pb (Front), size Medium, length 120 cm, in Green 17 and with MPure, MTex

Apron

	Front (Pb)	Length (Apron)	Article number* (Start)	(Colour code)	MPure (Inner material)	MFree (Inner material)
X-Small (XS)	0.25 mm Pb	90 cm	RA46020Z90	+CC	+55	+56
	0.25 mm Pb	100 cm	RA46020Z10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA46020Z11	+CC	+55	+56
	0.35 mm Pb	90 cm	RA46030Z90	+CC	+55	+56
	0.35 mm Pb	100 cm	RA46030Z10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA46030Z11	+CC	+55	+56
	0.50 mm Pb	90 cm	RA46050Z90	+CC	+55	+56
	0.50 mm Pb	100 cm	RA46050Z10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA46050Z11	+CC	+55	+56
Small (S)	0.25 mm Pb	90 cm	RA46020S90	+CC	+55	+56
	0.25 mm Pb	100 cm	RA46020S10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA46020S11	+CC	+55	+56
	0.35 mm Pb	90 cm	RA46030S90	+CC	+55	+56
	0.35 mm Pb	100 cm	RA46030S10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA46030S11	+CC	+55	+56
	0.50 mm Pb	90 cm	RA46050S90	+CC	+55	+56
	0.50 mm Pb	100 cm	RA46050S10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA46050S11	+CC	+55	+56
Medium (M)	0.25 mm Pb	100 cm	RA46020M10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA46020M11	+CC	+55	+56
	0.25 mm Pb	120 cm	RA46020M12	+CC	+55	+56
	0.35 mm Pb	100 cm	RA46030M10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA46030M11	+CC	+55	+56
	0.35 mm Pb	120 cm	RA46030M12	+CC	+55	+56
	0.50 mm Pb	100 cm	RA46050M10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA46050M11	+CC	+55	+56
	0.50 mm Pb	120 cm	RA46050M12	+CC	+55	+56
Large (L)	0.25 mm Pb	100 cm	RA46020L10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA46020L11	+CC	+55	+56
	0.25 mm Pb	120 cm	RA46020L12	+CC	+55	+56
	0.35 mm Pb	100 cm	RA46030L10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA46030L11	+CC	+55	+56
	0.35 mm Pb	120 cm	RA46030L12	+CC	+55	+56
	0.50 mm Pb	100 cm	RA46050L10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA46050L11	+CC	+55	+56
	0.50 mm Pb	120 cm	RA46050L12	+CC	+55	+56
X-Large (XL)	0.25 mm Pb	100 cm	RA46020X10	+CC	+55	+56
	0.25 mm Pb	110 cm	RA46020X11	+CC	+55	+56
	0.25 mm Pb	120 cm	RA46020X12	+CC	+55	+56
	0.35 mm Pb	100 cm	RA46030X10	+CC	+55	+56
	0.35 mm Pb	110 cm	RA46030X11	+CC	+55	+56
	0.35 mm Pb	120 cm	RA46030X12	+CC	+55	+56
	0.50 mm Pb	100 cm	RA46050X10	+CC	+55	+56
	0.50 mm Pb	110 cm	RA46050X11	+CC	+55	+56
	0.50 mm Pb	120 cm	RA46050X12	+CC	+55	+56



Colour codes



Inky Blue	51
New Red	52
Green 17	53
Modern Grey	54

Optional

Sewn-on shoulder border with embroidery (one line of text)
 Breast pocket (outer regular insertion pocket) without embroidery
 Breast pocket (outer regular insertion pocket) with embroidery
 Breast pocket (outer regular insertion pocket) with logo embroidery

Cap, MTex Outer material

Lead equivalent 0.25 mm Pb / 0.35 mm Pb / 0.50 mm Pb

Cap**Protection (Pb)****Size****Article number
(Start)****MPure
(Inner material)****MFree
(Inner material)**0.25 mm Pb
0.35 mm Pb
0.50 mm Pbuniversal
universal
universalRA410254
RA410354
RA410554+55
+55
+55+56
+56
+56**Thyroid protection, MTex** Outer material

Lead equivalent 0.50 mm Pb

Thyroid protection**Protection (Pb)****Closure****Article number
(Start)****MPure
(Inner material)****MFree
(Inner material)**

0.50 mm Pb

Hook&burr

RA414554KL

+55

+56

Bolero, MTex Outer material

Lead equivalent 0.25 mm Pb / 0.35 mm Pb / 0.50 mm Pb

Bolero**Protection (Pb)****Size****Article number
(Start)****MPure
(Inner material)****MFree
(Inner material)**

XS

0.25 mm Pb
0.35 mm Pb
0.50 mm PbX-Small
X-Small
X-SmallRA4372Z54
RA4373Z54
RA4375Z54+55
+55
+55+56
+56
+56

S

0.25 mm Pb
0.35 mm Pb
0.50 mm PbSmall
Small
SmallRA4372S54
RA4373S54
RA4375S54+55
+55
+55+56
+56
+56

M

0.25 mm Pb
0.35 mm Pb
0.50 mm PbMedium
Medium
MediumRA4372M54
RA4373M54
RA4375M54+55
+55
+55+56
+56
+56

L

0.25 mm Pb
0.35 mm Pb
0.50 mm PbLarge
Large
LargeRA4372L54
RA4373L54
RA4375L54+55
+55
+55+56
+56
+56

XL

0.25 mm Pb
0.35 mm Pb
0.50 mm PbX-Large
X-Large
X-LargeRA4372X54
RA4373X54
RA4375X54+55
+55
+55+56
+56
+56

Special hanger

Usage, capacity	Article number
For one item of clothing	AW106

3 Versions

for two, three or six items of clothing

Wall system

Usage, capacity	Article number
System for two special hangers	AW41502
System for three special hangers	AW41503
System for six special hangers	AW41506

Mobile system

for six or twelve items of clothing

Mobile system

Usage, capacity	Article number
For six items of clothing	AW42506
For twelve items of clothing	AW42512



MAVIG GmbH

Stahlgruberring 5
81829 Munich
Germany

Phone	+49 (0) 89 420 96 0
Fax	+49 (0) 89 420 96 200
e-Mail	info@mavig.com

MAVIG Nordic

Stockholm
Sweden

Phone	+46 (0) 722 25 25 68
e-Mail	larsson@mavig.com

MAVIG B.V.

Den Haag
The Netherlands

Phone	+31 (0) 70 33 11 688
Mobil	+31 (0) 61 595 43 48
e-Mail	simmonds@mavig.nl

MAVIG France SARL

Paris
France

Phone	+33 (0)1 30 59 46 23
e-Mail	info@mavig.fr

MAVIG Southeast Europe

Ljubljana
Slovenia

Phone	+386 (40) 6 33 900
e-Mail	dejak@mavig.com

Ti-Ba Enterprises, Inc.

Rochester, NY
USA

Phone	+1 (1) 585 247 1212
e-Mail	mavigusteam@mavig.com

Mena Medical Development

Beirut
Lebanon

Phone	+961 14 23 499
e-Mail	mavig@mena-md.com

AO Sante Medical Systems

Moscow
Russia

Phone	+7 499 551 55 73/75
e-Mail	info@sante.ru

